



Ex0

# Excel Formula<sup>®</sup> Orthodontic Denture Base Material

Our Orthodontic Material has been manufactured to a super-fine particle size which eliminates slumping during build-up. This material excels in the Salt-and-Pepper technique and with its ease of control gives technicians that consistent result they desire!

- ◆ Self-cure resin with superior clarity
- ◆ Super-fine particle size eliminates slumping during buildup
- ◆ 10 min. quick cure in hydroflask
- ◆ Powder mixes easily into liquid minimizing porosity
- ◆ Buffs to a high surface sheen within minutes
- ◆ High level of strength



## Available Products:



Lab Pack - 1 Kilo (2.2 lbs) powder, 500 mL (1 pint) clear liquid  
■ 12118 with clear liquid



1 Kilo Powder Only (2.2 lbs.)  
■ 12210 ■ 12230 LWT (Longer Working Time)



10 Kilo Powder Only (22 lbs.)  
■ 12220 ■ 12240 LWT (Longer Working Time)



1 Liter Liquid Only  
■ 12348 Clear



4x1 Liter Liquid Only  
■ 12358 Clear



Clear

